

Dry Pumps for Mass Spectrometry



Discover all the benefits of switching from a wet to a dry pump.



Zero contamination

FREE OF OIL SMELLS AND
PARTICULATES



Lower power

NXRI: 450W



Lower heat emission

NO AC UPGRADE REQUIRED



5 years interval

NO RESOURCE FOR
ANNUAL OIL CHANGES



High pumping density

ONE DRY PUMP CAN REPLACE
TWO WET PUMPS

You could reduce
your energy costs
up to 70% !



Wet pump replacement: two for one !

Liquid Chromatography Mass Spectrometry (LC-MS) integration example :

Before

2 X TYPICAL 40 M³/H



WET 2 DRY

After

1 X NXRI 90 M³/H



SCIEX
TRIPLE QUAD 6500



What about Gas Chromatography Mass Spectrometry (GC-MS) ?

PRODUCT PORTFOLIO

Our **nXDS range (scroll) pumps from Edwards** or the **VACUU·PURE from Vacuubrand** is ideal for this configuration.

CONTACT US TO CHOOSE THE RIGHT PUMP FOR YOUR SPECTROMETER

info@plasmadial.com | +41 32 926 26 06

